

ABX Minocal



LOT

CX 494

CAL



(Exp.) 2024-06-05
(YYYY - MM - DD)

Rev 1

| PARAMETRES PARAMETERS | UNITES UNITS | ABX Lysebio | | | | | | | | TOLERANCES TOLERANCE |
|--------------------------|---|----------------------------|----------------------------|---------|--|--|--|--|--|-------------------------|
| | | YUMIZEN | | YUMIZEN | | | | | | |
| | | H1500 H2500 HELO 1.x | H1500 H2500 HELO 2.0 | | | | | | | |
| GB WBC | 10 ⁹ /mm ³ ;10 ⁹ /l | 8.60 | 8.55 | | | | | | | ± 0.20 |
| GR RBC | 10 ⁶ /mm ³ ;10 ¹² /l | 4.52 | 4.53 | | | | | | | ± 0.06 |
| HB HGB | g/dl | 13.3 | 13.5 | | | | | | | ± 0.2 |
| | g/l | 133 | 135 | | | | | | | ± 2 |
| | mmol/l | 8.26 | 8.38 | | | | | | | ± 0.12 |
| PLA PLT | 10 ³ /mm ³ ;10 ⁹ /l | 243 | 243 | | | | | | | ± 10 |
| VGM MCV | µm ³ ; fl | 83.0 | 84.5 | | | | | | | ± 4.0 |

1300116192
Ref: TEMP-0387 Rev.51