



Yumizen G DDi 2 for YUMIZEN G1500/G1550



Reference No.:		1300036391	
Lot No.:		940224 (rev. 1)	
Expiry date:		2025-02-28	
Coagulometers		MASTER CURVE	
Optical	Yumizen G1500/G1550	µg(FEU)/mL	ΔOD
		4,00	0.2639
		3,00	0.2180
		2,00	0.1340
		1,00	0.0496
		0,50	0.0236
		0,25	0.0128
Slope value		0,06	

Yumizen G DDi 2 for YUMIZEN G800



Reference No.:		1300036391	
Lot No.:		940224 (rev. 1)	
Expiry date:		2025-02-28	
Coagulometers		MASTER CURVE	
Optical	Yumizen G800	µg(FEU)/mL	ΔOD
		4,00	0.2547
		3,00	0.1951
		2,00	0.1250
		1,00	0.0562
		0,50	0.0261
		0,25	0.0122
Slope value		0,06	



Reference No.:		1300036391	
Lot No.:		940224 (rev. 1)	
Expiry date:		2025-02-28	
Coagulometers		MASTER CURVE	
		µg(FEU)/mL	ΔOD
Optical	Yumizen G400 DDi/G200	4,00	0.1050
		3,00	0.0759
		2,00	0.0468
		1,00	0.0218
		0,50	0.0111
		0,25	0.0060
		Slope value	0,02

Important notice

For Yumizen G400 DDi and Yumizen G200, make sure that the end time is set to 180 seconds in the D-Dimer setup and the curve type **Lin-Lin p-p** is selected in **Calibration > D-Dimer > ug FEU/ml**.
 For more information, please contact your local HORIBA Medical representative.

Yumizen G DDi 2 for YUMIZEN G400 DDi/G200