

ABX Minocal



LOT

CX 494

CAL

 (Exp.) 2024-06-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x			YUMIZEN since V3						
		H550	H500 OT H500 CT		H500 OT	H500 CT H550					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.85	8.85		8.88	8.88					± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51	4.51		4.45	4.45					± 0.06
HB HGB	g/dl	13.0	13.0		13.0	13.0					± 0.2
	g/l	130	130		130	130					± 2
	mmol/l	8.07	8.07		8.07	8.07					± 0.12
HT HCT	%	37.9	37.9		36.8	36.6					± 1.0
	l/l	0.379	0.379		0.368	0.366					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	242	242		250	267					± 10
VMP MPV	µm ³ ; fl	10.3	10.3		N/A	N/A					± 0.5

1300113555
BACK / VERSO
Ref. TEMP-0387 Rev.51