

ABX Minocal



LOT		CX 491		CAL		(Exp.) 2024-03-05 (YYYY - MM - DD)							
		Rev 1		WHITEDIFF									
PARAMETRES PARAMETERS		UNITES UNITS		YUMIZEN V1.0 to V2.x		YUMIZEN since V3						TOLERANCES TOLERANCE	
				H550	H500 OT H500 CT		H500 OT	H500 CT H550					
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	8.75	8.75		8.70	8.70						± 0.20
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.57	4.57		4.57	4.57						± 0.06
HB	HGB	g/dl	12.8	12.8		12.9	12.9						± 0.2
		g/l	128	128		129	129						± 2
		mmol/l	7.95	7.95		8.01	8.01						± 0.12
HT	HCT	%	37.5	37.5		36.7	36.3						± 1.0
		l/l	0.375	0.375		0.367	0.363						± 0.010
PLA	PLT	10 ³ /mm ³ ;10 ⁹ /l	248	248		261	268						± 10
VMP	MPV	µm ³ ; fl	10.6	10.6		N/A	N/A						± 0.5

1300113555
BACK / VERSO
Ref. TEMP-0387 Rev.51