

ABX Minocal



LOT

CX 489

CAL

 (Exp.) 2024-01-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x			YUMIZEN since V3						
		H550	H500 OT H500 CT		H500 OT	H500 CT H550					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.80	8.80		8.80	8.80					± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.50		4.45	4.45					± 0.06
HB HGB	g/dl	13.0	13.0		13.1	13.1					± 0.2
	g/l	130	130		131	131					± 2
	mmol/l	8.07	8.07		8.14	8.14					± 0.12
HT HCT	%	37.8	37.8		36.5	36.4					± 1.0
	l/l	0.378	0.378		0.365	0.364					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	247	247		252	259					± 10
VMP MPV	µm ³ ; fl	10.7	10.7		N/A	N/A					± 0.5

1300113555
BACK / VERSO
Ref. TEMP-0387 Rev.50