

ABX Minocal



LOT

CX 473

Rev 1

CAL



(Exp.) 2022-09-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE	
		YUMIZEN V1.0 to V2.x				YUMIZEN since V3					
		H550	H500 OT H500 CT			H500 OT	H500 CT H550				
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.00	9.00			9.04	9.04				± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51	4.51			4.42	4.42				± 0.06
HB HGB	g/dl	12.8	12.8			12.8	12.8				± 0.2
	g/l	128	128			128	128				± 2
	mmol/l	7.95	7.95			7.97	7.97				± 0.12
HT HCT	%	37.4	37.4			35.8	35.8				± 1.0
	l/l	0.374	0.374			0.358	0.358				± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	260	260			269	269				± 10
VMP MPV	μm ³ ; fl	10.3	10.3			10.4	10.4				± 0.5

Ref: TEMP-0387 Rev.48 BACK / VERSO 1300105208