

LOT		CX 483		CAL		(Exp.) 2023-07-05 (YYYY - MM - DD)						
PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE	
		YUMIZEN V1.0 to V2.x		YUMIZEN since V3								
		H550	H500 OT H500 CT		H500 OT	H500 CT H550						
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.90	8.90		8.90	8.90					± 0.20	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.52		4.50	4.50					± 0.06	
HB HGB	g/dl	12.8	12.8		12.8	12.8					± 0.2	
	g/l	128	128		128	128					± 2	
	mmol/l	7.95	7.95		7.95	7.95					± 0.12	
HT HCT	%	37.1	37.1		35.5	34.8					± 1.0	
	l/l	0.371	0.371		0.355	0.348					± 0.010	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	258	258		283	288					± 10	
VMP MPV	μm ³ ; fl	10.3	10.3		N/A	N/A					± 0.5	

Ref: TEMP-0387 Rev.49

BACK / VERSO

1300113555

 1300113555
 BACK / VERSO
 Ref. TEMP-0387 Rev.49