

ABX Minocal



(Exp.) 2023-05-05

(YYYY - MM - DD)

LOT		CX 481		CAL		(Exp.) 2023-05-05 (YYYY - MM - DD)									
PARAMETRES PARAMETERS		UNITES UNITS		WHITEDIFF								TOLERANCES TOLERANCE			
				YUMIZEN V1.0 to V2.x				YUMIZEN since V3							
				H550	H500 OT H500 CT		H500 OT	H500 CT H550							
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	8.80	8.80		8.80	8.80					± 0.20			
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42	4.42		4.45	4.45					± 0.06			
HB	HGB	g/dl	12.5	12.5		12.5	12.5					± 0.2			
		g/l	125	125		125	125					± 2			
		mmol/l	7.76	7.76		7.76	7.76					± 0.12			
HT	HCT	%	36.5	36.5		35.2	34.4					± 1.0			
		l/l	0.365	0.365		0.352	0.344					± 0.010			
PLA	PLT	10 ³ /mm ³ ;10 ⁹ /l	248	248		275	278					± 10			
VMP	MPV	μm ³ ; fl	10.7	10.7		9.2	9.2					± 0.5			

 1300105208
 BACK / VERSO
 Ref. TEMP-0387 Rev.48