

ABX Minocal

LOT		CX 479		CAL		(Exp.) 2023-03-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x			YUMIZEN since V3						
		H550	H500 OT H500 CT		H500 OT	H500 CT H550					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.80	8.80		8.75	8.80					± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.60	4.60		4.55	4.48					± 0.06
HB HGB	g/dl	13.0	13.0		12.9	12.9					± 0.2
	g/l	130	130		129	129					± 2
	mmol/l	8.07	8.07		8.01	8.01					± 0.12
HT HCT	%	38.2	38.2		36.5	35.3					± 1.0
	l/l	0.382	0.382		0.365	0.353					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	259	259		257	285					± 10
VMP MPV	µm ³ ; fl	11.6	11.6		9.2	11.2					± 0.5

Ref: TEMP-0387 Rev.48

1300105208

BACK / VERSO

 1300105208
 BACK / VERSO
 Ref. TEMP-0387 Rev.48