

ABX Minocal


LOT
CX 478
CAL

(Exp.) 2023-02-05

(YYYY - MM - DD)

 Rev **2**

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x			YUMIZEN since V3						
		H550	H500 OT H500 CT		H500 OT	H500 CT H550					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.80	8.80		8.75	8.75					± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.59	4.59		4.50	4.50					± 0.06
HB HGB	g/dl	13.1	13.1		13.0	13.3					± 0.2
	g/l	131	131		130	133					± 2
	mmol/l	8.14	8.14		8.07	8.26					± 0.12
HT HCT	%	38.3	38.3		35.3	35.7					± 1.0
	l/l	0.383	0.383		0.353	0.357					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	249	249		264	274					± 10
VMP MPV	μm ³ ; fl	11.0	11.0		9.2	9.2					± 0.5

 1300105208
 BACK / VERSO
 Ref. TEMP-0387 Rev.48