

ABX Minocal



LOT

CX 478

CAL



(Exp.) 2023-02-05

Rev **2**

(YYYY - MM - DD)

| PARAMETRES PARAMETERS | UNITES UNITS | WHITEDIFF | | | | | | | | TOLERANCES TOLERANCE |
|--------------------------|---|----------------------|--------------------|--|--|------------------|-----------------|--|--|-------------------------|
| | | YUMIZEN V1.0 to V2.x | | | | YUMIZEN since V3 | | | | |
| | | H550 | H500 OT H500 CT | | | H500 OT | H500 CT H550 | | | |
| GB WBC | 10 ³ /mm ³ ;10 ⁹ /l | 8.80 | 8.80 | | | 8.75 | 8.75 | | | ± 0.20 |
| GR RBC | 10 ⁶ /mm ³ ;10 ¹² /l | 4.59 | 4.59 | | | 4.50 | 4.50 | | | ± 0.06 |
| HB HGB | g/dl | 13.1 | 13.1 | | | 13.0 | 13.3 | | | ± 0.2 |
| | g/l | 131 | 131 | | | 130 | 133 | | | ± 2 |
| | mmol/l | 8.14 | 8.14 | | | 8.07 | 8.26 | | | ± 0.12 |
| HT HCT | % | 38.3 | 38.3 | | | 35.3 | 35.7 | | | ± 1.0 |
| | l/l | 0.383 | 0.383 | | | 0.353 | 0.357 | | | ± 0.010 |
| PLA PLT | 10 ³ /mm ³ ;10 ⁹ /l | 249 | 249 | | | 264 | 274 | | | ± 10 |
| VMP MPV | µm ³ ; fl | 11.0 | 11.0 | | | 9.2 | 9.2 | | | ± 0.5 |

Ref: TEMP-0387 Rev.48 BACK / VERSO 1300105208