

ABX Minocal

LOT		CX 476		CAL		(Exp.) 2022-12-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x		YUMIZEN since V3							
		H550	H500 OT H500 CT		H500 OT	H500 CT H550					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.93	8.93		8.56	9.26					± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.54		4.45	4.57					± 0.06
HB HGB	g/dl	12.9	12.9		12.6	13.0					± 0.2
	g/l	129	129		126	130					± 2
	mmol/l	8.01	8.01		7.82	8.07					± 0.12
HT HCT	%	37.7	37.7		35.0	35.8					± 1.0
	l/l	0.377	0.377		0.350	0.358					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	249	249		259	265					± 10
VMP MPV	μm ³ ; fl	11.5	11.5		10.6	10.8					± 0.5

 1300105208
 BACK / VERSO
 Ref. TEMP-0387 Rev.48