

ABX Minocal



LOT

CX 495

CAL



(Exp.) 2024-07-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x				YUMIZEN since V3				
		H550	H500 OT H500 CT			H500 OT	H500 CT H550			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.90	8.90			8.90	8.90			± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.50			4.47	4.47			± 0.06
HB HGB	g/dl	12.9	12.9			13.0	13.0			± 0.2
	g/l	129	129			130	130			± 2
	mmol/l	8.01	8.01			8.07	8.07			± 0.12
HT HCT	%	37.4	37.4			36.6	35.2			± 1.0
	l/l	0.374	0.374			0.366	0.352			± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	242	242			244	257			± 10
VMP MPV	μm ³ ; fl	9.8	9.8			N/A	N/A			± 0.5

Ref: TEMP-0387 Rev.51 BACK / VERSO 1300113555