

ABX Minocal



LOT		CX 490		CAL		(Exp.) 2024-02-05 (YYYY - MM - DD)							
		Rev 1		WHITEDIFF									
PARAMETRES PARAMETERS		UNITES UNITS		YUMIZEN V1.0 to V2.x		YUMIZEN since V3						TOLERANCES TOLERANCE	
				H550	H500 OT H500 CT		H500 OT	H500 CT H550					
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	8.85	8.85		8.85	8.85						± 0.20
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.47	4.47		4.47	4.47						± 0.06
HB	HGB	g/dl	12.6	12.6		12.8	12.8						± 0.2
		g/l	126	126		128	128						± 2
		mmol/l	7.82	7.82		7.95	7.95						± 0.12
HT	HCT	%	36.9	36.9		34.8	35.0						± 1.0
		l/l	0.369	0.369		0.348	0.350						± 0.010
PLA	PLT	10 ³ /mm ³ ;10 ⁹ /l	252	252		263	272						± 10
VMP	MPV	µm ³ ; fl	10.2	10.2		N/A	N/A						± 0.5

1300113555
BACK / VERSO
Ref. TEMP-0387 Rev.51