

ABX Minocal

LOT		CX 487		CAL		(Exp.) 2023-11-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x			YUMIZEN since V3						
		H550	H500 OT H500 CT		H500 OT	H500 CT H550					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.10	9.10		9.10	9.10					± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.61	4.61		4.52	4.52					± 0.06
HB HGB	g/dl	13.0	13.0		13.0	13.0					± 0.2
	g/l	130	130		130	130					± 2
	mmol/l	8.07	8.07		8.07	8.07					± 0.12
HT HCT	%	38.0	38.0		35.9	35.2					± 1.0
	l/l	0.380	0.380		0.359	0.352					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	242	242		264	269					± 10
VMP MPV	μm ³ ; fl	10.3	10.3		N/A	N/A					± 0.5

Ref: TEMP-0387 Rev.50

BACK / VERSO

1300113555

 1300113555
 BACK / VERSO
 Ref. TEMP-0387 Rev.50