

LOT		CX 484		CAL		(Exp.) 2023-08-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x			YUMIZEN since V3						
		H550	H500 OT H500 CT		H500 OT	H500 CT H550					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.80	8.80		8.80	8.80					± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.57	4.57		4.54	4.54					± 0.06
HB HGB	g/dl	13.2	13.2		13.3	13.3					± 0.2
	g/l	132	132		133	133					± 2
	mmol/l	8.20	8.20		8.26	8.26					± 0.12
HT HCT	%	38.6	38.6		36.9	36.4					± 1.0
	l/l	0.386	0.386		0.369	0.364					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	253	253		272	277					± 10
VMP MPV	μm ³ ; fl	10.4	10.4		N/A	N/A					± 0.5

Ref: TEMP-0387 Rev.49

1300113555

BACK / VERSO

 1300113555
 BACK / VERSO
 Ref. TEMP-0387 Rev.49