

LOT		CX 482		CAL		(Exp.) 2023-06-05 (YYYY - MM - DD)									
PARAMETRES PARAMETERS		UNITES UNITS		WHITEDIFF								TOLERANCES TOLERANCE			
				YUMIZEN V1.0 to V2.x				YUMIZEN since V3							
				H550	H500 OT H500 CT		H500 OT	H500 CT H550							
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	9.00	9.00		9.00	9.00					± 0.20			
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42	4.42		4.43	4.43					± 0.06			
HB	HGB	g/dl	12.8	12.8		12.8	12.8					± 0.2			
		g/l	128	128		128	128					± 2			
		mmol/l	7.95	7.95		7.95	7.95					± 0.12			
HT	HCT	%	37.1	37.1		36.1	35.1					± 1.0			
		l/l	0.371	0.371		0.361	0.351					± 0.010			
PLA	PLT	10 ³ /mm ³ ;10 ⁹ /l	255	255		275	280					± 10			
VMP	MPV	µm ³ ; fl	10.4	10.4		9.2	9.2					± 0.5			

Ref: TEMP-0387 Rev.48

1300105208

BACK / VERSO

 1300105208
 BACK / VERSO
 Ref. TEMP-0387 Rev.48