

ABX Minocal



LOT		CX 495		CAL		(Exp.) 2024-07-05 (YYYY - MM - DD)							
		Rev 1		WHITEDIFF									
PARAMETRES PARAMETERS	UNITES UNITS	YUMIZEN V1.0 to V2.x			YUMIZEN since V3						TOLERANCES TOLERANCE		
		H550	H500 OT H500 CT		H500 OT	H500 CT H550							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.90	8.90		8.90	8.90					± 0.20		
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.50		4.47	4.47					± 0.06		
HB HGB	g/dl	12.9	12.9		13.0	13.0					± 0.2		
	g/l	129	129		130	130					± 2		
	mmol/l	8.01	8.01		8.07	8.07					± 0.12		
HT HCT	%	37.4	37.4		36.6	35.2					± 1.0		
	l/l	0.374	0.374		0.366	0.352					± 0.010		
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	242	242		244	257					± 10		
VMP MPV	µm ³ ; fl	9.8	9.8		N/A	N/A					± 0.5		

1300113555
BACK / VERSO
Ref. TEMP-0387 Rev.51