

# ABX Minocal

| LOT                      |                                                       | CX 472               |                    | CAL |                  | <div><div></div><div>(Exp.) 2022-08-05</div><div>( YYYY - MM - DD )</div></div> |  |  |  |  |                         |
|--------------------------|-------------------------------------------------------|----------------------|--------------------|-----|------------------|---------------------------------------------------------------------------------|--|--|--|--|-------------------------|
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                       | WHITEDIFF            |                    |     |                  |                                                                                 |  |  |  |  | TOLERANCES<br>TOLERANCE |
|                          |                                                       | YUMIZEN V1.0 to V2.x |                    |     | YUMIZEN since V3 |                                                                                 |  |  |  |  |                         |
|                          |                                                       | H550                 | H500 OT<br>H500 CT |     | H500 OT          | H500 CT<br>H550                                                                 |  |  |  |  |                         |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 8.90                 | 8.90               |     | 9.00             | 9.00                                                                            |  |  |  |  | ± 0.20                  |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | 4.57                 | 4.57               |     | 4.48             | 4.48                                                                            |  |  |  |  | ± 0.06                  |
| HB HGB                   | g/dl                                                  | 13.2                 | 13.2               |     | 13.3             | 13.3                                                                            |  |  |  |  | ± 0.2                   |
|                          | g/l                                                   | 132                  | 132                |     | 133              | 133                                                                             |  |  |  |  | ± 2                     |
|                          | mmol/l                                                | 8.20                 | 8.20               |     | 8.26             | 8.26                                                                            |  |  |  |  | ± 0.12                  |
| HT HCT                   | %                                                     | 38.6                 | 38.6               |     | 35.0             | 35.0                                                                            |  |  |  |  | ± 1.0                   |
|                          | l/l                                                   | 0.386                | 0.386              |     | 0.350            | 0.350                                                                           |  |  |  |  | ± 0.010                 |
| PLA PLT                  | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 256                  | 256                |     | 255              | 255                                                                             |  |  |  |  | ± 10                    |
| VMP MPV                  | µm <sup>3</sup> ; fl                                  | 10.9                 | 10.9               |     | 10.6             | 10.6                                                                            |  |  |  |  | ± 0.5                   |

Ref: TEMP-0387 Rev.48 BACK / VERSO 1300105208