

# ABX Minocal

| LOT                      |                                                       | CX 485               |                    | CAL |                  | <div><div></div><div>(Exp.) 2023-09-05</div><div>(YYYY - MM - DD)</div></div> |  |  |  |  |                         |
|--------------------------|-------------------------------------------------------|----------------------|--------------------|-----|------------------|-------------------------------------------------------------------------------|--|--|--|--|-------------------------|
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                       | WHITEDIFF            |                    |     |                  |                                                                               |  |  |  |  | TOLERANCES<br>TOLERANCE |
|                          |                                                       | YUMIZEN V1.0 to V2.x |                    |     | YUMIZEN since V3 |                                                                               |  |  |  |  |                         |
|                          |                                                       | H550                 | H500 OT<br>H500 CT |     | H500 OT          | H500 CT<br>H550                                                               |  |  |  |  |                         |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 9.10                 | 9.10               |     | 9.10             | 9.10                                                                          |  |  |  |  | ± 0.20                  |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | 4.58                 | 4.58               |     | 4.53             | 4.53                                                                          |  |  |  |  | ± 0.06                  |
| HB HGB                   | g/dl                                                  | 13.0                 | 13.0               |     | 13.0             | 13.0                                                                          |  |  |  |  | ± 0.2                   |
|                          | g/l                                                   | 130                  | 130                |     | 130              | 130                                                                           |  |  |  |  | ± 2                     |
|                          | mmol/l                                                | 8.07                 | 8.07               |     | 8.07             | 8.07                                                                          |  |  |  |  | ± 0.12                  |
| HT HCT                   | %                                                     | 38.0                 | 38.0               |     | 36.6             | 36.1                                                                          |  |  |  |  | ± 1.0                   |
|                          | l/l                                                   | 0.380                | 0.380              |     | 0.366            | 0.361                                                                         |  |  |  |  | ± 0.010                 |
| PLA PLT                  | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 246                  | 246                |     | 275              | 270                                                                           |  |  |  |  | ± 10                    |
| VMP MPV                  | μm <sup>3</sup> ; fl                                  | 10.0                 | 10.0               |     | N/A              | N/A                                                                           |  |  |  |  | ± 0.5                   |

Ref: TEMP-0387 Rev.49

1300113555

BACK / VERSO

 1300113555  
 BACK / VERSO  
 Ref. TEMP-0387 Rev.49