

ABX Minocal

LOT		CX 488		CAL		(Exp.) 2023-12-05 (YYYY - MM - DD)						
PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE	
		YUMIZEN V1.0 to V2.x				YUMIZEN since V3						
		H550	H500 OT H500 CT		H500 OT	H500 CT H550						
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.95	8.95		8.98	8.98					± 0.20	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.41	4.41		4.39	4.39					± 0.06	
HB HGB	g/dl	12.7	12.7		12.7	12.7					± 0.2	
	g/l	127	127		127	127					± 2	
	mmol/l	7.89	7.89		7.89	7.89					± 0.12	
HT HCT	%	37.0	37.0		35.8	35.2					± 1.0	
	l/l	0.370	0.370		0.358	0.352					± 0.010	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	239	239		250	261					± 10	
VMP MPV	μm ³ ; fl	10.2	10.2		N/A	N/A					± 0.5	

Ref: TEMP-0387 Rev.50

BACK / VERSO

1300113555

 1300113555
 BACK / VERSO
 Ref. TEMP-0387 Rev.50