

ABX Minocal

LOT		CX 476		CAL		(Exp.) 2022-12-05 (YYYY - MM - DD)									
PARAMETRES PARAMETERS		UNITES UNITS		WHITEDIFF								TOLERANCES TOLERANCE			
				YUMIZEN V1.0 to V2.x				YUMIZEN since V3							
				H550	H500 OT H500 CT		H500 OT	H500 CT H550							
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	8.93	8.93		8.56	9.26					± 0.20			
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.54		4.45	4.57					± 0.06			
HB	HGB	g/dl	12.9	12.9		12.6	13.0					± 0.2			
		g/l	129	129		126	130					± 2			
		mmol/l	8.01	8.01		7.82	8.07					± 0.12			
HT	HCT	%	37.7	37.7		35.0	35.8					± 1.0			
		l/l	0.377	0.377		0.350	0.358					± 0.010			
PLA	PLT	10 ³ /mm ³ ;10 ⁹ /l	249	249		259	265					± 10			
VMP	MPV	μm ³ ; fl	11.5	11.5		10.6	10.8					± 0.5			

Ref: TEMP-0387 Rev.48

BACK / VERSO

1300105208

 1300105208
 BACK / VERSO
 Ref. TEMP-0387 Rev.48