

ABX Minocal

LOT
CX 474
CAL

(Exp.) 2022-10-05

(YYYY - MM - DD)

 Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x			YUMIZEN since V3						
		H550	H500 OT H500 CT		H500 OT	H500 CT H550					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.85	8.85		8.74	9.07					± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.54		4.45	4.58					± 0.06
HB HGB	g/dl	12.9	12.9		12.8	13.0					± 0.2
	g/l	129	129		128	130					± 2
	mmol/l	8.01	8.01		7.95	8.07					± 0.12
HT HCT	%	37.7	37.7		36.2	36.6					± 1.0
	l/l	0.377	0.377		0.362	0.366					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	253	253		254	279					± 10
VMP MPV	μm ³ ; fl	10.4	10.4		10.7	10.8					± 0.5

 1300105208
 BACK / VERSO
 Ref. TEMP-0387 Rev.48