

ABX Minocal



LOT		CX 495		CAL		(Exp.) 2024-07-05 (YYYY - MM - DD)							
		Rev 1		WHITEDIFF									
				YUMIZEN V1.0 to V2.x				YUMIZEN since V3					
PARAMETRES PARAMETERS		UNITES UNITS		H550	H500 OT H500 CT		H500 OT	H500 CT H550					TOLERANCES TOLERANCE
GB WBC		10 ³ /mm ³ ;10 ⁹ /l		8.90	8.90		8.90	8.90					± 0.20
GR RBC		10 ⁶ /mm ³ ;10 ¹² /l		4.50	4.50		4.47	4.47					± 0.06
HB HGB		g/dl		12.9	12.9		13.0	13.0					± 0.2
		g/l		129	129		130	130					± 2
		mmol/l		8.01	8.01		8.07	8.07					± 0.12
HT HCT		%		37.4	37.4		36.6	35.2					± 1.0
		l/l		0.374	0.374		0.366	0.352					± 0.010
PLA PLT		10 ³ /mm ³ ;10 ⁹ /l		242	242		244	257					± 10
VMP MPV		µm ³ ; fl		9.8	9.8		N/A	N/A					± 0.5

1300113555
BACK / VERSO
Ref. TEMP-0387 Rev.51