

ABX Minocal



(Exp.) 2023-09-05

(YYYY - MM - DD)

LOT		CX 485		CAL		(Exp.) 2023-09-05 (YYYY - MM - DD)									
PARAMETRES PARAMETERS		UNITES UNITS		WHITEDIFF								TOLERANCES TOLERANCE			
				YUMIZEN V1.0 to V2.x				YUMIZEN since V3							
				H550	H500 OT H500 CT		H500 OT	H500 CT H550							
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	9.10	9.10		9.10	9.10					± 0.20			
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.58	4.58		4.53	4.53					± 0.06			
HB	HGB	g/dl	13.0	13.0		13.0	13.0					± 0.2			
		g/l	130	130		130	130					± 2			
		mmol/l	8.07	8.07		8.07	8.07					± 0.12			
HT	HCT	%	38.0	38.0		36.6	36.1					± 1.0			
		l/l	0.380	0.380		0.366	0.361					± 0.010			
PLA	PLT	10 ³ /mm ³ ;10 ⁹ /l	246	246		275	270					± 10			
VMP	MPV	μm ³ ; fl	10.0	10.0		N/A	N/A					± 0.5			

 1300113555
 BACK / VERSO
 Ref. TEMP-0387 Rev.49