

ABX Minocal



LOT		CX 491		CAL		(Exp.) 2024-03-05 (YYYY - MM - DD)							
		Rev 1		WHITEDIFF									
PARAMETRES PARAMETERS	UNITES UNITS	YUMIZEN V1.0 to V2.x			YUMIZEN since V3						TOLERANCES TOLERANCE		
		H550	H500 OT H500 CT		H500 OT	H500 CT H550							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.75	8.75		8.70	8.70					± 0.20		
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.57	4.57		4.57	4.57					± 0.06		
HB HGB	g/dl	12.8	12.8		12.9	12.9					± 0.2		
	g/l	128	128		129	129					± 2		
	mmol/l	7.95	7.95		8.01	8.01					± 0.12		
HT HCT	%	37.5	37.5		36.7	36.3					± 1.0		
	l/l	0.375	0.375		0.367	0.363					± 0.010		
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	248	248		261	268					± 10		
VMP MPV	μm ³ ; fl	10.6	10.6		N/A	N/A					± 0.5		

1300113555
BACK / VERSO
Ref. TEMP-0387 Rev.51