

ABX Minocal



(Exp.) 2023-12-05

(YYYY - MM - DD)

LOT		CX 488		CAL		(Exp.) 2023-12-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x			YUMIZEN since V3						
		H550	H500 OT H500 CT		H500 OT	H500 CT H550					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.95	8.95		8.98	8.98					± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.41	4.41		4.39	4.39					± 0.06
HB HGB	g/dl	12.7	12.7		12.7	12.7					± 0.2
	g/l	127	127		127	127					± 2
	mmol/l	7.89	7.89		7.89	7.89					± 0.12
HT HCT	%	37.0	37.0		35.8	35.2					± 1.0
	l/l	0.370	0.370		0.358	0.352					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	239	239		250	261					± 10
VMP MPV	μm ³ ; fl	10.2	10.2		N/A	N/A					± 0.5

 1300113555
 BACK / VERSO
 Ref. TEMP-0387 Rev.50