

ABX Minocal

LOT		CX 477		CAL		(Exp.) 2023-01-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x			YUMIZEN since V3						
		H550	H500 OT H500 CT		H500 OT	H500 CT H550					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.78	8.78		8.45	9.16					± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.62	4.62		4.54	4.66					± 0.06
HB HGB	g/dl	13.6	13.6		13.4	13.7					± 0.2
	g/l	136	136		134	137					± 2
	mmol/l	8.45	8.45		8.32	8.51					± 0.12
HT HCT	%	39.3	39.3		37.3	38.3					± 1.0
	l/l	0.393	0.393		0.373	0.383					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	270	270		255	275					± 10
VMP MPV	µm ³ ; fl	11.1	11.1		10.6	11.0					± 0.5

Ref: TEMP-0387 Rev.48

BACK / VERSO

1300105208

 1300105208
 BACK / VERSO
 Ref. TEMP-0387 Rev.48