

ABX Minocal



LOT

CX 486

CAL



(Exp.) 2023-10-05

Rev **2**

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x				YUMIZEN since V3				
		H550	H500 OT H500 CT			H500 OT	H500 CT H550			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.80	8.80			8.80	8.80			± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.52			4.50	4.50			± 0.06
HB HGB	g/dl	12.9	12.9			12.9	12.9			± 0.2
	g/l	129	129			129	129			± 2
	mmol/l	8.01	8.01			8.01	8.01			± 0.12
HT HCT	%	37.7	37.7			36.2	35.7			± 1.0
	l/l	0.377	0.377			0.362	0.357			± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	243	243			263	265			± 10
VMP MPV	μm ³ ; fl	9.7	9.7			N/A	N/A			± 0.5

Ref: TEMP-0387 Rev.50 BACK / VERSO 1300113555