

ABX Minocal

LOT		CX 484		CAL		(Exp.) 2023-08-05 (YYYY - MM - DD)							
PARAMETRES PARAMETERS		UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE		
			YUMIZEN V1.0 to V2.x				YUMIZEN since V3						
			H550	H500 OT H500 CT		H500 OT	H500 CT H550						
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	8.80	8.80		8.80	8.80					± 0.20	
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.57	4.57		4.54	4.54					± 0.06	
HB	HGB	g/dl	13.2	13.2		13.3	13.3					± 0.2	
		g/l	132	132		133	133					± 2	
		mmol/l	8.20	8.20		8.26	8.26					± 0.12	
HT	HCT	%	38.6	38.6		36.9	36.4					± 1.0	
		l/l	0.386	0.386		0.369	0.364					± 0.010	
PLA	PLT	10 ³ /mm ³ ;10 ⁹ /l	253	253		272	277					± 10	
VMP	MPV	μm ³ ; fl	10.4	10.4		N/A	N/A					± 0.5	

Ref: TEMP-0387 Rev.49

1300113555

BACK / VERSO

 1300113555
 BACK / VERSO
 Ref. TEMP-0387 Rev.49