

ABX Minocal


LOT
CX 482
CAL

(Exp.) 2023-06-05

(YYYY - MM - DD)

 Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x			YUMIZEN since V3						
		H550	H500 OT H500 CT		H500 OT	H500 CT H550					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.00	9.00		9.00	9.00					± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42	4.42		4.43	4.43					± 0.06
HB HGB	g/dl	12.8	12.8		12.8	12.8					± 0.2
	g/l	128	128		128	128					± 2
	mmol/l	7.95	7.95		7.95	7.95					± 0.12
HT HCT	%	37.1	37.1		36.1	35.1					± 1.0
	l/l	0.371	0.371		0.361	0.351					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	255	255		275	280					± 10
VMP MPV	μm ³ ; fl	10.4	10.4		9.2	9.2					± 0.5

 1300105208
 BACK / VERSO
 Ref. TEMP-0387 Rev.48