

ABX Minocal


LOT
CX 481

 Rev **2**
CAL

(Exp.) 2023-05-05

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x			YUMIZEN since V3						
		H550	H500 OT H500 CT		H500 OT	H500 CT H550					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.80	8.80		8.80	8.80					± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42	4.42		4.45	4.45					± 0.06
HB HGB	g/dl	12.5	12.5		12.5	12.5					± 0.2
	g/l	125	125		125	125					± 2
	mmol/l	7.76	7.76		7.76	7.76					± 0.12
HT HCT	%	36.5	36.5		35.2	34.4					± 1.0
	l/l	0.365	0.365		0.352	0.344					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	248	248		275	278					± 10
VMP MPV	μm ³ ; fl	10.7	10.7		9.2	9.2					± 0.5

 1300105208
 BACK / VERSO
 Ref. TEMP-0387 Rev.48