

ABX Minocal



LOT

CX 479

CAL



(Exp.) 2023-03-05

Rev **2**

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x				YUMIZEN since V3				
		H550	H500 OT H500 CT			H500 OT	H500 CT H550			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.80	8.80			8.75	8.80			± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.60	4.60			4.55	4.48			± 0.06
HB HGB	g/dl	13.0	13.0			12.9	12.9			± 0.2
	g/l	130	130			129	129			± 2
	mmol/l	8.07	8.07			8.01	8.01			± 0.12
HT HCT	%	38.2	38.2			36.5	35.3			± 1.0
	l/l	0.382	0.382			0.365	0.353			± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	259	259			257	285			± 10
VMP MPV	µm ³ ; fl	11.6	11.6			9.2	11.2			± 0.5

Ref: TEMP-0387 Rev.48 BACK / VERSO 1300105208