

ABX Minocal

LOT		CX 474		CAL		(Exp.) 2022-10-05 (YYYY - MM - DD)						
PARAMETRES PARAMETERS		UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE	
			YUMIZEN V1.0 to V2.x		YUMIZEN since V3							
			H550	H500 OT H500 CT		H500 OT	H500 CT H550					
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	8.85	8.85		8.74	9.07					± 0.20
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.54		4.45	4.58					± 0.06
HB	HGB	g/dl	12.9	12.9		12.8	13.0					± 0.2
		g/l	129	129		128	130					± 2
		mmol/l	8.01	8.01		7.95	8.07					± 0.12
HT	HCT	%	37.7	37.7		36.2	36.6					± 1.0
		l/l	0.377	0.377		0.362	0.366					± 0.010
PLA	PLT	10 ³ /mm ³ ;10 ⁹ /l	253	253		254	279					± 10
VMP	MPV	μm ³ ; fl	10.4	10.4		10.7	10.8					± 0.5

Ref: TEMP-0387 Rev.48

1300105208

BACK / VERSO

 1300105208
 BACK / VERSO
 Ref. TEMP-0387 Rev.48