

ABX Minocal

LOT		CX 486		CAL		(Exp.) 2023-10-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x			YUMIZEN since V3						
		H550	H500 OT H500 CT		H500 OT	H500 CT H550					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.80	8.80		8.80	8.80					± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.52		4.50	4.50					± 0.06
HB HGB	g/dl	12.9	12.9		12.9	12.9					± 0.2
	g/l	129	129		129	129					± 2
	mmol/l	8.01	8.01		8.01	8.01					± 0.12
HT HCT	%	37.7	37.7		36.2	35.7					± 1.0
	l/l	0.377	0.377		0.362	0.357					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	243	243		263	265					± 10
VMP MPV	µm ³ ; fl	9.7	9.7		N/A	N/A					± 0.5

Ref: TEMP-0387 Rev.50

BACK / VERSO

1300113555

 1300113555
 BACK / VERSO
 Ref. TEMP-0387 Rev.50