

LOT		CX 484		CAL		(Exp.) 2023-08-05 (YYYY - MM - DD)									
PARAMETRES PARAMETERS		UNITES UNITS		WHITEDIFF								TOLERANCES TOLERANCE			
				YUMIZEN V1.0 to V2.x				YUMIZEN since V3							
				H550	H500 OT H500 CT		H500 OT	H500 CT H550							
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	8.80	8.80		8.80	8.80					± 0.20			
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.57	4.57		4.54	4.54					± 0.06			
HB	HGB	g/dl	13.2	13.2		13.3	13.3					± 0.2			
		g/l	132	132		133	133					± 2			
		mmol/l	8.20	8.20		8.26	8.26					± 0.12			
HT	HCT	%	38.6	38.6		36.9	36.4					± 1.0			
		l/l	0.386	0.386		0.369	0.364					± 0.010			
PLA	PLT	10 ³ /mm ³ ;10 ⁹ /l	253	253		272	277					± 10			
VMP	MPV	μm ³ ; fl	10.4	10.4		N/A	N/A					± 0.5			

Ref: TEMP-0387 Rev.49

BACK / VERSO

1300113555

 1300113555
 BACK / VERSO
 Ref. TEMP-0387 Rev.49