

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MC 432

CONTROL



(Exp.) 2024-10-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS		UNITES UNITS		ABX Lysebio													
				CONTROL			1	CONTROL			2	CONTROL			3	TOLERANCES TOLERANCE	
				MICROS		PENTRA	TOLERANCES TOLERANCE	MICROS		PENTRA	TOLERANCES TOLERANCE	MICROS		PENTRA			
				EMI CRP	LC-767G	MS CRP		EMI CRP	LC-767G	MS CRP		EMI CRP	LC-767G	MS CRP			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	2.0	2.0	2.1	± 0.4	7.7		7.6	7.8	± 1.0		20.4	20.4	21.0	± 2.2		
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	2.38	2.37	2.33	± 0.15	4.49	4.46	4.49	± 0.18	5.50	5.47	5.54	± 0.25				
HB HGB	g/dl	6.3	6.3	6.0	± 0.4	13.5	13.3	13.4	± 0.5	17.5	17.4	17.7	± 0.7				
	g/l	63	63	60	± 4	135	133	134	± 5	175	174	177	± 7				
	mmol/l	3.91	3.91	3.73	± 0.25	8.38	8.26	8.32	± 0.31	10.87	10.81	10.99	± 0.44				
HT HCT	%	17.1	17.1	16.5	± 2.0	37.3	36.6	36.6	± 2.5	48.4	47.6	47.9	± 3.0				
	l/l	0.171	0.171	0.165	± 0.020	0.373	0.366	0.366	± 0.025	0.484	0.476	0.479	± 0.030				
VGM MCV	µm ³ ; fl	72	72	71.0	± 4	83	82	81.5	± 4	88	87	86.5	± 4				
TGMH MCH	pg	26.5	26.6	25.8	± 2.4	30.1	29.8	29.8	± 2.4	31.8	31.8	31.9	± 2.8				
	fmol	1.64	1.65	1.60	± 0.15	1.87	1.85	1.85	± 0.15	1.98	1.98	1.98	± 0.18				
CCMH MCHC	g/dl	36.8	36.9	36.3	± 3.0	36.2	36.4	36.6	± 3.0	36.2	36.6	36.9	± 3.0				
	g/l	368	369	363	± 30	362	364	366	± 30	362	366	369	± 30				
	mmol/l	22.83	22.93	22.52	± 1.86	22.50	22.58	22.74	± 1.86	22.45	22.71	22.94	± 1.86				
IDR RDW	%	14.5	14.0	14.5	± 3.0	14.0	13.0	13.0	± 3.0	13.5	12.5	11.5	± 3.0				
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	66	73	63	± 20	245	243	240	± 40	468	477	507	± 60				
VPM MPV	µm ³ ; fl	7.3	7.9	8.7	± 2.0	7.9	8.0	8.6	± 2.0	7.3	7.4	7.7	± 2.0				
LYMPHO	#	1.3	1.3	N/A	± 0.2	2.6	2.5	N/A	± 0.5	3.5	3.4	N/A	± 1.1				
	%	66.0	65.5	N/A	± 8.0	33.5	33.0	N/A	± 6.0	17.0	16.5	N/A	± 5.0				
MONO	#	0.2	0.2	N/A	± 0.2	0.5	0.5	N/A	± 0.5	1.0	1.1	N/A	± 1.0				
	%	8.0	8.0	N/A	± 7.2	6.5	7.0	N/A	± 6.3	5.0	5.5	N/A	± 5.0				
GRANULO	#	0.5	0.5	N/A	± 0.2	4.6	4.6	N/A	± 0.8	15.9	15.9	N/A	± 2.0				
	%	26.0	26.5	N/A	± 10.0	60.0	60.0	N/A	± 10.0	78.0	78.0	N/A	± 10.0				
CRP	mg/dl	0.49	0.44	0.44	± 0.40	2.78	2.86	2.64	± 0.80	5.01	4.87	5.10	± 1.50				
	mg/l	4.9	4.4	4.4	± 4.0	27.8	28.6	26.4	± 8.0	50.1	48.7	51.0	± 15.0				

Ref: TEMP-1070 Rev.43 BACKVERSO 1300092170