

ABX Minotrol Retic



LOT RX 448

CONTROL

(Exp.) **2024-09-05**
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL 1			TOLERANCES TOLERANCE	CONTROL 2			TOLERANCES TOLERANCE	CONTROL 3			TOLERANCES TOLERANCE
		YUMIZEN				YUMIZEN				YUMIZEN			
		H2500 HELO 2.0				H2500 HELO 2.0				H2500 HELO 2.0			
Méthode Manuelle / Manual													
RET	%	1.3			± 0.7	5.6			± 2.3	10.6			± 3.2
RET	%	1.72			± 0.86	5.02			± 1.26	9.83			± 1.97
RET	10 ⁶ /mm ³	0.060			± 0.030	0.173			± 0.043	0.344			± 0.069
	10 ⁹ /l	60			± 30	173			± 43	344			± 69
CRC	%	1.13			± 0.57	3.26			± 0.98	6.18			± 1.55
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.47			± 0.30	3.45			± 0.30	3.50			± 0.30
IRF	--	0.22			± 0.13	0.26			± 0.16	0.28			± 0.17
RHCC	pg	21.6			± 4.3	15.6			± 3.1	10.8			± 2.7

Ref: TEMP-0382 Rev.45 BACK / VERSO 1300116193