

ABX Minocal



LOT		CX 496		CAL		(Exp.) 2024-08-05 (YYYY - MM - DD)						
		Rev 1		ABX Lysebio								
PARAMETRES PARAMETERS	UNITES UNITS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60		60 60 C+ ES60	XL80		
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.3	9.1	9.4	9.2	9.5	9.2		9.2	9.2	± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42	4.38	4.37	4.40	4.46	4.43		4.41	4.40	± 0.06	
HB HGB	g/dl	13.0	12.9	13.0	13.1	13.1	13.0		13.1	13.1	± 0.2	
	g/l	130	129	130	131	131	130		131	131	± 2	
	mmol/l	8.07	8.01	8.07	8.14	8.14	8.07		8.14	8.14	± 0.12	
HT HCT	%	36.9	36.4	35.2	35.9	35.9	35.4		35.9	35.9	± 1.0	
	l/l	0.369	0.364	0.352	0.359	0.359	0.354		0.359	0.359	± 0.010	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	239	236	227	233	242	233		236	233	± 10	
VMP MPV	µm ³ ; fl	8.5	8.7	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5	

1300090308
BACK / VERSO
Ref. TEMP-0387 Rev.51