

ABX Minocal



LOT		CX 496		CAL		 (Exp.) 2024-08-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN								
		H1500	H1500								
		H2500	H2500								
		HELO 1.x	HELO 2.0								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.70	8.82								± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.44	4.49								± 0.06
HB HGB	g/dl	13.0	13.3								± 0.2
	g/l	130	133								± 2
	mmol/l	8.07	8.26								± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	243	239								± 10
VGM MCV	µm ³ ; fl	82.5	83.8								± 4.0

1300116192
Ref. TEMP-0387 Rev.51