



HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée
 B.P. 7290
 34184 Montpellier cedex 4
 FRANCE
 www.horiba-abx.com



A03A00009B

Yumizen G PT Liq 4 for Yumizen G Line



D 4 9 3 0 5 0 9 3 6 3 7 3 1 1 2 4 D



D 5 1 0 0 1 1 9 5 0 0 1 6 7 3 3 0 D



D 6 2 1 5 2 5 0 2 6 1 1 1 9 1 2 0 D



D 7 1 5 3 D

Reference No.:		1300036373		
Lot No.:		930509 (rev. 1)		
Expiry date:		2024-11-30		
Coagulometers		ISI	MASTER CURVE	
Optical	Yumizen G Line		1,20	%
		100		11,9
		50		16,7
		33		21,5
		25		26,1
	MNPT	11,9		
Important notice				
<p>The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.</p>				