

ABX Minocal

LOT

CX 497

CAL



(Exp.) 2024-09-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0				9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.49				4.46					± 0.06
HB HGB	g/dl	13.1				13.1					± 0.2
	g/l	131				131					± 2
	mmol/l	8.14				8.14					± 0.12
HT HCT	%	37.0				37.9					± 1.0
	l/l	0.370				0.379					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	223				233					± 10
VMP MPV	μm ³ ; fl	9.3				9.7					± 0.5

ABX Minocal

LOT

CX 497

CAL



(Exp.) 2024-09-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS	MICROS							
			Scil Vet ABC Plus Scil Vet ABC ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		8.7							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.45							± 0.06
HB HGB	g/dl		12.8							± 0.2
	g/l		128							± 2
	mmol/l		7.95							± 0.12
HT HCT	%		36.9							± 1.0
	l/l		0.369							± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l		233							± 10
VMP MPV *	µm ³ ; fl		9.1							± 0.5

Ref: TEMP-0387 Rev.52 FRONT / RECTO 1300113555

* N/A for Micros Scil Vet ABC