

ABX Minocal



LOT		CX 497		CAL				(Exp.) 2024-09-05 (YYYY - MM - DD)				
		Rev 1		ABX Lysebio								
PARAMETRES PARAMETERS	UNITES UNITS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60		60 60 C+ ES60	XL80		
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.0	9.3	9.1	9.3	9.1		9.2	9.1	± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.47	4.44	4.44	4.47	4.51	4.45		4.50	4.47	± 0.06	
HB HGB	g/dl	13.0	12.9	13.1	13.2	13.1	13.0		13.1	13.2	± 0.2	
	g/l	130	129	131	132	131	130		131	132	± 2	
	mmol/l	8.07	8.01	8.14	8.20	8.14	8.07		8.14	8.20	± 0.12	
HT HCT	%	37.1	36.4	35.7	36.0	36.1	35.6		36.2	36.0	± 1.0	
	l/l	0.371	0.364	0.357	0.360	0.361	0.356		0.362	0.360	± 0.010	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	234	239	229	227	237	226		234	227	± 10	
VMP MPV	µm ³ ; fl	9.0	9.1	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5	

1300090308
BACK / VERSO
Ref: TEMP-0387 Rev.52