

ABX Minocal



LOT

CX 497

CAL



(Exp.) 2024-09-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN						
		H1500	H1500	H2500	H2500					
		HELO 1.x	HELO 2.0							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.60	8.60							± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.54							± 0.06
HB HGB	g/dl	13.0	13.0							± 0.2
	g/l	130	130							± 2
	mmol/l	8.07	8.07							± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	241	241							± 10
VGM MCV	µm ³ ; fl	81.9	81.9							± 4.0

1300116192
Ref: TEMP-0387 Rev.52