

ABX Minocal



LOT		CX 497		CAL		(Exp.) 2024-09-05 (YYYY - MM - DD)					
Rev 1		ABX Lysebio									
PARAMETRES PARAMETERS	UNITES UNITS	YUMIZEN	YUMIZEN								TOLERANCES TOLERANCE
		H1500	H1500								
		H2500	H2500								
		HELO 1.x	HELO 2.0								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.60	8.60								± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.54								± 0.06
HB HGB	g/dl	13.0	13.0								± 0.2
	g/l	130	130								± 2
	mmol/l	8.07	8.07								± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	241	241								± 10
VGM MCV	µm ³ ; fl	81.9	81.9								± 4.0

1300116192
Ref. TEMP-0387 Rev.52