

ABX Minocal



LOT		CX 498		CAL		(Exp.) 2024-10-05 (YYYY - MM - DD)					
Rev 1		ABX Lysebio									
PARAMETRES PARAMETERS	UNITES UNITS	YUMIZEN	YUMIZEN								TOLERANCES TOLERANCE
		H1500	H1500								
		H2500	H2500								
		HELO 1.x	HELO 2.0								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.50	8.50								± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.48	4.48								± 0.06
HB HGB	g/dl	12.8	12.8								± 0.2
	g/l	128	128								± 2
	mmol/l	7.95	7.95								± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	230	230								± 10
VGM MCV	µm ³ ; fl	81.3	81.3								± 4.0

1300116192
Ref. TEMP-0387 Rev.52