

ABX Minocal



LOT

CX 499

CAL



(Exp.) 2024-11-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60 Care ST				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9				8.9					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.39				4.42					± 0.06
HB HGB	g/dl	12.7				12.7					± 0.2
	g/l	127				127					± 2
	mmol/l	7.89				7.89					± 0.12
HT HCT	%	36.4				36.2					± 1.0
	l/l	0.364				0.362					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	236				234					± 10
VMP MPV	µm ³ ; fl	8.9				7.9					± 0.5

Ref: TEMP-0387 Rev.52 FRONT / RECTO 1300090308