

# ABX Minocal



LOT

CX 499

CAL

 (Exp.) 2024-11-05  
(YYYY - MM - DD)

Rev 1

| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                       | ABX Lysebio |         |        |        |                      |        |        |                     |        | TOLERANCES<br>TOLERANCE |
|--------------------------|-------------------------------------------------------|-------------|---------|--------|--------|----------------------|--------|--------|---------------------|--------|-------------------------|
|                          |                                                       | MICROS      | MICROS  | PENTRA | PENTRA | PENTRA               | PENTRA | PENTRA | PENTRA              | PENTRA |                         |
|                          |                                                       | EMI CRP     | LC-767G | MS CRP | XLR    | DX NEXUS<br>DF NEXUS | MS60   |        | 60<br>60 C+<br>ES60 | XL80   |                         |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 9.0         | 8.9     | 9.1    | 8.8    | 9.0                  | 8.8    |        | 8.8                 | 8.8    | ± 0.2                   |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | 4.40        | 4.40    | 4.37   | 4.43   | 4.43                 | 4.40   |        | 4.43                | 4.43   | ± 0.06                  |
| HB HGB                   | g/dl                                                  | 12.5        | 12.7    | 12.6   | 12.7   | 12.5                 | 12.5   |        | 12.8                | 12.7   | ± 0.2                   |
|                          | g/l                                                   | 125         | 127     | 126    | 127    | 125                  | 125    |        | 128                 | 127    | ± 2                     |
|                          | mmol/l                                                | 7.76        | 7.89    | 7.82   | 7.89   | 7.76                 | 7.76   |        | 7.95                | 7.89   | ± 0.12                  |
| HT HCT                   | %                                                     | 35.6        | 35.9    | 34.3   | 35.0   | 34.3                 | 34.1   |        | 35.2                | 35.0   | ± 1.0                   |
|                          | l/l                                                   | 0.356       | 0.359   | 0.343  | 0.350  | 0.343                | 0.341  |        | 0.352               | 0.350  | ± 0.010                 |
| PLA PLT                  | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 234         | 238     | 229    | 226    | 234                  | 230    |        | 232                 | 226    | ± 10                    |
| VMP MPV                  | µm <sup>3</sup> ; fl                                  | 8.8         | 8.9     | N/A    | N/A    | N/A                  | N/A    |        | N/A                 | N/A    | ± 0.5                   |

1300090308  
BACK / VERSO  
Ref: TEMP-0387 Rev.52