

ABX Minocal



LOT		CX 499		CAL		(Exp.) 2024-11-05 (YYYY - MM - DD)						
		Rev 1										
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE	
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA		
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60		60 60 C+ ES60	XL80		
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	8.9	9.1	8.8	9.0	8.8		8.8	8.8	± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.40	4.40	4.37	4.43	4.43	4.40		4.43	4.43	± 0.06	
HB HGB	g/dl	12.5	12.7	12.6	12.7	12.5	12.5		12.8	12.7	± 0.2	
	g/l	125	127	126	127	125	125		128	127	± 2	
	mmol/l	7.76	7.89	7.82	7.89	7.76	7.76		7.95	7.89	± 0.12	
HT HCT	%	35.6	35.9	34.3	35.0	34.3	34.1		35.2	35.0	± 1.0	
	l/l	0.356	0.359	0.343	0.350	0.343	0.341		0.352	0.350	± 0.010	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	234	238	229	226	234	230		232	226	± 10	
VMP MPV	µm ³ ; fl	8.8	8.9	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5	

1300090308
BACK / VERSO
Ref: TEMP-0387 Rev.52