

ABX Minocal



LOT

CX 499

CAL



(Exp.) 2024-11-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN						
		H1500	H1500	H2500	H2500					
		HELO 1.x	HELO 2.0							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.40	8.40							± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.50							± 0.06
HB HGB	g/dl	12.6	12.6							± 0.2
	g/l	126	126							± 2
	mmol/l	7.82	7.82							± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	235	235							± 10
VGM MCV	µm ³ ; fl	79.1	79.1							± 4.0

1300116192
Ref: TEMP-0387 Rev.52